

# TOTOKU

1.3 Megapixel  
18.1" Color Display

# CCL182 /r

## CCL Series

## 1.3MP color LCD with color calibration function

Dual analog/digital interfaces with a color calibration function.  
Its built-in luminance stabilizing system provides stability in luminance.



*Flat Display System  
for Medical Imaging*

# 1.3

Megapixel  
Color Display

- 18.1"
- Analog & Digital
- High luminance 270cd/m<sup>2</sup>
- Contrast 400:1
- Calibration function
- 10Bit LUT
- Multi function Stand

