

TOTOKU

2 Megapixel
20.1" Color Display

CCL204

CCL Series

2MP color LCD with color calibration function

Dual analog/digital interfaces with a color calibration function.
Its built-in luminance stabilizing system provides stability in luminance.

*Flat Display System
for Medical Imaging*



2 Megapixel
Color Display

20.1"	Analog & Digital	High luminance 300cd/m ²	Contrast 1000:1	Calibration function	10Bit LUT	Multi function Stand
-------	------------------	-------------------------------------	-----------------	----------------------	-----------	----------------------

