

TOTOKU


2 Megapixel
21.3" Color Display

CCL212 /r

CCL Series

2MP color LCD with color calibration function and high luminance of 450cd/m²

Dual analog/digital interfaces with a color calibration function.
Its built-in luminance stabilizing system provides stability in luminance.



Flat Display System
for Medical Imaging

2

Megapixel
Color Display

High Luminance

21.3"

Analog
& Digital

High
luminance
450cd/m²

Contrast
450:1

Calibration
function

10Bit
LUT

LED
Indicator

Multi
function
Stand

