

TOTOKU

2 Megapixel
20.1" Monochrome Display


ME201L /r

ME Series

2MP monochrome LCD with calibration function and high contrast ratio of 1000:1

Dual analog/digital interface model best suited for various imaging modalities and the PACS. The built-in calibration function enables DICOM GSDF compliant grayscale presentation.

Flat Display System
for Medical Imaging



2

Megapixel
Monochrome Display

- 20.1"
- Analog & Digital
- High luminance 700cd/m²
- Contrast 1000:1
- Calibration function
- 10Bit LUT
- Multi function Stand

