

TOTOKU

2 Megapixel
21.3" Monochrome Display

ME251i


ME251i / C (clear base), ME251i / B (blue base)

High-luminance, long life model with easy evaluation of DICOM GSDF conformance.

Combined with a viewer, simultaneous display of 1024 shades of gray in compliant with DICOM GSDF is possible. The color/monochrome conversion function converts color images into crisp and clear monochrome images.



Higher Image Quality
& Total Management



2

Megapixel
Monochrome Display



21.3"

Digital Input

High luminance 1800cd/m²

Contrast 700:1

Calibration function

10-bit display

Color/Monochrome Conversion

OSD

LED Indicator

Protective filter

Multi function Stand

